



FRASTEMA

ophthalmic equipments

88MC OPTOMETRIST DELUXE



General

Height	1305mm
Width	580mm
Depth	880mm
Rotation	45 deg. right / 45 deg. Left

Head Rest

Vertical Range	Fixed
----------------	-------

Seat

Horizontal Range	65mm
Vertical Range	180mm
Minimum Height	525mm
Width	460mm
Depth	480mm

Back Rest

Backward	60 deg.
Height	880mm
Width	470mm

Foot Rest

Back Rotation	90 deg.
---------------	---------

Arm Rests

Back Rotation	90 deg.
---------------	---------

Weight	70 kg
---------------	-------

Distributed by


"Focussed on Eyecare"

Hanson Instruments Ltd,
44 Washford Industrial Estate, Heming Road,
Redditch, Worcestershire B98 0DP

01527 501077

www.hansoninstruments.co.uk