

hanson TW-1002 Manual Focimeter



- Internal reading Focimeter
- Classic manual focimeter with internal reading design.
- Green LED light source with large visual field.
- Installed with Prism compensator.
- Ease of use whether right and left hander
- Ergonomic design optimized operational feel

Sphere lens	-25.00D~+25.00D
Step	0.25D(-5D~+5D,Step 0.125D)
Axis	0°~180°,Step1°
Prism basal angel	0°~180°,Step1°
Prism degree	Standard Type0~5Δ,Step0.5Δ
	Enhancement Type 0~20Δ,Step1Δ
Lens diameter	ø16mm~ø80mm
Eyepiece visibility adjusting range	-5.00D~+5.00D
Size and weight	280(W)×150(D)×450(H)mm/4.9kg
Power supply	Input AC 100-240V, 50/60Hz ;Output DC 6V 2A