

TW-1004 Manual Focimeter



- Green LED light source
- Ranges of measurement: 0 to $\pm 25D$ Minimum: 0.125D graduation up to $\pm 5D$ 0.25D graduations over $\pm 5D$ to $\pm 25D$
- Cylinder Axis Range: 0 \sim 180° step 1°
- With Prism Compensator
- Prismatic Power Range: Prism diopter of Pattern A: at 1 \square interval in the range of 0 \sim 5 \square Prism diopter of Pattern B: at 1 \square interval in the range of 0 \sim 20 \square
- Prism base angle :at 50 intervals in the range of 0 \sim 180 degree
- Eyepiece Adjustment: $\pm 5D$
- Measurand Lens Diameter: $\phi 16 \sim \phi 80mm$
- Dimensions: 275mm(h) \times 130mm(w) \times 455mm(d)
- Weight: 5.6kg
- Lamp: 220V15W 110V15W