


# Polaris

## Tear-Film Analyzer



# With Polaris you can...

Display and assess tear film lipid patterns



Evaluate tear meniscus

Fit a contact lens and evaluate the tear lipid pattern after fitting

Measure Tear break up time

The basic principle is to analyze the tear-film patterns generated by the interference between the white light and the tear-film itself (fringes).

Depending on the thickness and composition of the tearfilm, fringes will show different shapes and colors.



44 Washford Industrial Estate, Heming Road

Redditch. B98 0DP, UK

Tel : 0044 (0)1527 501077