



1-800-979-2020 • 610-889-0200 • FAX 610-889-3233 • [info@accutome.com](mailto:info@accutome.com) • [www.accutome.com](http://www.accutome.com)

## Total Solutions Provider for all of your ophthalmic needs.

### ACCUTOME SURGICAL

Diamond Knives, Steel Knives, Sutures and Handheld Surgical Instruments.

### ACCUTOME CLINICAL

Office Medications and Supplies.

### ACCUTOME DIAGNOSTIC

A-Scans, B-Scans, Pachymeters, Tonometer and UBM.

### ACCUTOME REPAIR CENTER

- Quick turnaround for most brands of Diamond Knives and Handheld Instruments.
- Diagnostic Equipment.

### ORDERING

Call toll free (800) 979.2020 US and Canada  
(610) 889.0200 International  
+31 (0) 485.350300 Europe  
Monday thru Friday 8:30AM - 8PM EST  
Fax: (610) 889.3233  
Email: [info@accutome.com](mailto:info@accutome.com)  
Web: [www.accutome.com](http://www.accutome.com)

### SHIPPING

Ground - Via UPS Ground Service  
Air - Via best method

Saturday - Orders must be placed before 3PM EST  
Friday to ensure Saturday delivery.

Payment should be made by check, Bank Transfer, Visa,  
MasterCard or American Express.

## CONTACT US



### WATCH RESOURCES

[youtube.com/c/accutome](https://youtube.com/c/accutome)



### TWEET US

[twitter.com/accutome](https://twitter.com/accutome)



### GOOGLE+

[plus.google.com/+Accutome/](https://plus.google.com/+Accutome/)



### LIKE US

[facebook.com/pages/Accutome-Inc](https://facebook.com/pages/Accutome-Inc)



### FOLLOW US

[linkedin.com/company/accutome-inc](https://linkedin.com/company/accutome-inc)



Distributed by:



## Tonometer & Pachymeters

**ACCUPEN®** HANDHELD TONOMETER

**PACHPEN®** HANDHELD PACHYETER

**ACCUPACH VI®** DESKTOP PACHYETER





## When the Pressure's On... The AccuPen® Delivers

**The AccuPen® - Handheld Applanation Tonometer**  
Product Number: 24-3000

### Accurate

- Gravity Offset Technology provides precise IOP measurements with less calibration compared to other handheld tonometers.
- Micro strain gauge technology combined with proprietary algorithm provides reproducible results.

### Easy-to-Use

- Ergonomic design combined with sighting lines allow for easy visualization of the cornea, facilitating both centralization and perpendicularity.
- Contoured molding with large, visible LCD allows the AccuPen® to be used ambidextrously and fit comfortably in any size hand.

### Ready-to-Use

- Eliminates the constant need for "up, down and around" calibration.
- Long lasting lithium battery lasts 4-6 months.



### AccuTip Covers

Sleeved, sanitized and individually packed  
**AX9950** - Box of 100  
**AX11932** - Box of 300

### Portability

- Weighs only 3.0 oz.
- Comes with custom, heavy duty traveling case and lanyard for safe transportation.

### Versatility

- Fast measurements with patient in any position.
- Frees up valuable office space.

### Smarter Software

- Stores and averages nine measurements.
- Built in IOP Correction Calculator adjusts IOP based on corneal thickness.

## The Next Generation in Pachymetry

**The PachPen® - Handheld Pachymeter**  
Product Number: 24-5100

Accutome's PachPen® is the perfect combination of accuracy, portability and comfort. The superior ergonomic design fits comfortably in any hand, offering a major improvement from other handheld models.

### Accurate

- 65 MHz sampling probe ensures precise measurements.
- Uses Digital Signal Analysis to analyze wave patterns and eliminate non perpendicular measurements.

### Easy-to-Use

- Ergonomic design, combined with sighting lines, allow for easy visualization of the cornea, facilitating both centralization and perpendicularity.
- Contoured molding with large, visible LCD allows the PachPen® to be used ambidextrously and fit comfortably in any size hand.

### Ready-to-Use

- Digital ultrasound analysis eliminates need for calibration.
- Long lasting lithium battery lasts 4-6 months.

### Portability

- Weighs only 3.0 oz.
- Comes with custom, heavy duty traveling case and lanyard for safe transportation.

### Versatility

- Fast measurements with patient in any position.
- Frees up valuable office space.

### Smarter Software

- Stores and averages nine measurements.
- Built in IOP Correction Calculator adjusts IOP based on corneal thickness.

## Specifications

### AccuPen® Specifications

Accuracy	+/- 2mmHg (compared to Goldman Tonometer)
Weight	3.0 oz.
Dimensions	7 ¼" x 1 ¼"
Transducer Dimensions	1.0 mm diameter
Measurement Range	5-60mmHg
Power Source	3.6 V Lithium Battery

### PachPen® and Pach VI® Specifications

Accuracy	+/- 5 Microns
Weight	3.0 oz.
Dimensions	7 ¼" x 1 ¼"
Probe Dimension	2.5 mm diameter
Measurement Range	300-999 microns
Power Source	3.6 V Lithium Battery (PachPen) Medical grade Power Supply (AccuPach VI)

## AccuPach VI®

**Desktop Pachymeter**  
Product Number: 24-6200

- **Easy-to-Use** - User-friendly touch-tone screen allows the user to master the AccuPach VI® in minutes.
- **Patients are #1** - Revolutionary voice output feature calls out readings, allowing the user to concentrate more on the cornea.
- **Superior Accuracy** - State-of-the-art digital technology (Accutome's Digital Signal Analysis), combined with a probe sampling of 65 MHz, ensures precise readings.
- **USB Interface** - Link to computer or printer via USB Memory Stick.
- **Confidence** - Digital waveform analysis helps to ensure measurements in seconds.
- **Flexibility** - Adjustable handle/stand allows viewing on different angles.

