



INTRODUCING THE...
JVL/1 LED
OPHTHALMOMETER



INTRODUCING THE...

JVL/1 LED

OPHTHALMOMETER

High quality performance
ophthalmometer for the examination of:

- Refractive capacity of the cornea
- Bending radius of the cornea
- Corneal astigmatism
- Axle direction of the two measured meridians

Wide surface corneal measuring (3,4mm).
Supplied with red and green filters which
provide a white colour on superimposed
parts, enabling precision and exact
measurements. The new joystick allows
the user to check all instrument
movements at the same time.



Standard accessories:

- Protective cover
- Control sphere with support for contact lenses 50D
- Chin rest paper
- Two protection fuses

Technical specification:

- Measuring range:
dioptries from 30,00 to 60,00 D -
radius from 5,6 to 11,3mm
- Resolution:
dioptries 0,25 D - radius 0,05mm
- Useful measuring diameter:
3,4 mm equal to 9 square mm
- Rating scale:
TABE & INTERNATIONAL
- Lines of sight:
according to Javal
- Instrument base:
cross-wise movement to be stopped
- Chin rest:
adjustable height with shutter device
- Feeding voltage:
120Vac/60Hz - 230 Vac/50Hz
- Instrument voltage: 6Vac
- LED illumination
- Fuses: 5x20mm, 120V - 250mA, 230V - 125mA
- Absorbing power (VA): 24VA

hanson

"Focussed on Eyecare"

**For more information,
please contact Hanson
Instruments, the
exclusive UK distributor
of CSO on 01527 501077**

44 Washford Industrial Estate,
Heming Road, Redditch B98 0DP
Web: www.hansoninstruments.co.uk